



## INSTRUCTIONS

This kit is intended for the collection of bovine fecal samples from feedlot pens to submit for Q•resist testing.

1. Select a pen and record the pen identification in the space provided on the included submission form.
2. Also record the Collection Date in the space provided (DD/MM/YYYY format).
3. Select a relatively fresh excrement sample (a sample recently produced by an animal is best) and using the spoon provided in the kit, collect a heaping spoonful and deposit the fecal material into one of the containers provided.
4. Repeat 'Step 3' collection from 4 more relatively fresh excrement samples in the pen.
5. Ensure lids have been closed on each sample container.
6. Complete the remainder of the submission form and send samples along with the submission form to:

Quantum Genetix  
101-110 Research Drive  
Saskatoon, SK S7N 3R3  
Canada

## Submission Form

1. Submitted by  (Please print clearly)

Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ Prov/State: \_\_\_\_\_ Postal/Zip Code: \_\_\_\_\_

Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_

Email: \_\_\_\_\_

Send results to (if different than above):

Name: \_\_\_\_\_

Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_

Email: \_\_\_\_\_

### 2. Sample Information

Pen ID # \_\_\_\_\_

Collection Date DD / MM / YYYY

Lab ID (lab use only)